Tedequim

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TEDEQUIM® REPARADOR BIOCERÁMICO

PÁGINA..... 2

PÁGINA..... 3

INSTRUCTIONS FOR USE

TEDEQUIM® BIOCERAMIC REPAIR

PAGE..... 4

PAGE..... 5

TEDEQUIM®

REPARADOR BIOCERÁMICO

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TEDEQUIM®

BIOCERAMIC REPAIR

BASED ON CALCIUM SILICATE TECHNOLOGY

Description

Bioceramic Repair Material (Putty) It is bioactive and based on tricalcium silicate technology. It promotes the repair and regeneration of dental tissues and provides an airtight seal. It has properties similar to those of healthy dentin and can replace dentin at the coronary and root levels.

Instructions for Use

- 1) Isolate the tooth to treat.
- 2) Take the bottle of powder and shake it before opening.
- 3) Remove one level scoop of powder and place it on a glass slab. Close the bottle immediately.
- 4) Shake the dropper bottle containing the liquid. Open it and place a drop near the powder.
- 5) Spatulate for about 30 seconds until a putty mixture forms.
- 6) Place the mixture on the area to be treated.
- 7) Immediately wash the materials used.
- 8) The mixture will harden in a few minutes.

Indication

- Treatment of root canal perforations.
- Sealing of internal and external communicating resorptions.
- Retro-obturation in paraendodontic surgeries.
- Direct pulp capping.
- Pulpotomy.
- Apexification.
- Apexogenesis.

Features

- Easy to use and insert into the cavity.
- Release of calcium ions induces tissue regeneration.
- It has bactericidal and bacteriostatic properties due to the high pH (approximately 12) generated by the hydration of tricalcium silicate.
- Good radiopacity due to zirconium oxide content.
- Its good expansion allows it to adhere to dentin and prevent bacterial filtration.
- It does not contain resin.
- It does not stain.
- It is not affected by a wet environment.
- The average particle size is less than 25 microns.

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Limit of application

- Restoration of important dentin losses subjected to strong pressures.
- Sthetic restoration of the anterior sector.
- Treatment of teeth with irreversible pulpitis.

Adverse reactions

There are currently no known adverse reactions.

Composition

- Powder: White and mainly composed of tricalcium silicate, dicalcium silicate, silicon oxide, calcium carbonate, and zirconium oxide.
- Liquid: Calcium chloride, a hydrophilic plasticizing agent, and water.

Contenido

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Preservation and storage

Store tightly closed, in the metallized bag together with the moisture absorber at a temperature not exceeding 30°C.

Autorized by ANMAT PM 1277-36

Technical director: Miryam Della Vedova - Pharmacist - Licence 2107

Exclusive sale to professionals and sanitary facilities

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