**Tedequim**

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**Clarident TA**  
**Micro-abrasive Technique**  
**Read all instructions before use**  
**EXCLUSIVE PROFESSIONAL USE**

**Clarident TA** is a complete kit for the application of the micro-abrasive technique to eliminate stains produced by fluorosis in enamel of moderate or severe degree, where the whitening with peroxides is not enough. In the most severe cases it allows to reduce the tone of the color considerably. This technique uses a paste based on silicon carbide and hydrochloric acid and can be combined with home bleaching treatment with carbamide peroxide.

A treatment with micro-abrasive causes a loss of enamel of approximately 0.2 microns, reason why it is considered more convenient than restorative techniques (veneers) that require whitening with peroxides -to lower the tone of the stain- prior to the application of the veneer.

Micro-abrasive is also indicated in situations such as hypoplastic lesions and irregularities in the texture of the enamel, for white spots resulting from inactive caries, paralyzed pigmented caries and other structural defects of the superficial enamel that are normally pigmented.

**Instructions for use:**

Make an accurate diagnosis of the origin of the abnormal coloration. It is not necessary to anesthetize.

1. Check that the gum is healthy in the area to be bleached.
2. Treat cavities and possible leaks. In case of sensitivity, apply desensitizing agent NaFRESIN or KDESIN before starting this treatment.
3. Remove bacterial plaque. It is convenient to carry out a complete prophylaxis treatment. Do not use materials that contains glycerin or fluorine.
4. Cover the soft tissue and lips with TEDEQUIM GINGIVAL AND LABIAL PROTECTOR.
5. Isolate those teeth to whiten with rubber dam and leave at least one tooth visible to control the coloration.
6. Prepare the abrasive paste to be used according to the instructions indicated below.
  - 6.1. Apply the paste with the help of the spatula on the tooth to be bleached and leave for 5 seconds.
  - 6.2. Rub with the spatula from top to bottom and from mesial to distal, following the location of the stain(s) only, for about 20 seconds or more. It can also be rubbed mechanically. It is not convenient to apply more than 60 seconds at a time.
  - 6.3. Remove the product by sucking or with a cotton swab and rinse with plenty of water and observe.
  - 6.4. Repeat point 6.1, 6.2 and 6.3 no more than 10 times until the desired coloration is achieved.
7. To ensure that the acid was completely eliminated after rinsing, wet a cotton swab with TEDEQUIM Neutralizing Solution and apply to the treated teeth for a few seconds. Then rinse with plenty of water.
8. Carry out a smooth polishing of the treated elements.
9. Finally use fluoride at neutral pH NEUTRODENT or KDESIN.

Notes: Whatever the instrument used in point 6.2 (manual or mechanical technique), the results obtained are similar in terms of micro-abrasive effectiveness, amount of removed enamel and alteration of the surface roughness of the tooth enamel. The difference between the techniques consists in the greater or less time necessary for the removal of the stains, leaving for the dentist the decision of the technique to be used.

If at the end of the procedure the professional considers necessary, a dental whitening treatment with peroxides, whether ambulatory or professional, can be carried out.

#### **Preparation of the abrasive paste:**

- a) Place a spoonful (level) in a Dappen glass or similar well-cleaned container.
- b) Add 2 drops of hydrochloric acid, mixing until obtaining a homogeneous paste. This amount depends on the measured powder, you can add up to 3 drops.

NOTE: The dose indicated above is reaches to one element, repeating the procedure up to 8 times.

Following the powder/liquid ratio, a greater quantity can be prepared according to the number of elements to be treated.

The complete kit reaches about 70 teeth.

#### **Recommendations:**

- Use gloves.
- Protect the patient's eyes, assistant and professional with protective glasses.
- After preparing the abrasive paste, keep the water and the suction pump handy in case of splashing.
- Do not leave the Dappen cup with the paste above the patient's face.
- Keep the product away from children.

#### **Caution:**

Prevent contact with skin or gingival tissue. In case this happens, neutralize and wash with large amount of water, for 15 minutes.

#### **Contraindications:**

Do not use in patients with tooth sensitivity. In that case, treat it previously with NaFRESIN or KDESIN before starting the abrasion technique.

Do not apply on dentin or decayed teeth. Only enamel.

Do not place in areas with possible leaks.

#### **Presentation:**

The kit contains:

- Gingival protector: 9 g.
- Hydrochloric acid 18%: 6 ml.
- Abrasive Powder: 15 g.
- Accessories.

**Storage:**

Keep between 4°C and 40 °C protected from light.

**Authorized by ANMAT PM 1277-24**

**Technical director: Miryam Della Vedova - Pharmacist - License 2107**

**Exclusive sale to professionals and sanitary facilities**

**Manufactured by TEDEQUIM SRL**

**Bv. De los Polacos 6136, Los Boulevares.**

**X5147GGP Córdoba, Argentina.**

**Tel: + 54 3543 448260**

**www.tedequim.com**

**Rev. C**