Tedequim

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www.tedequim.com

Rev. B

CLARIDENT X EXCLUSIVE PROFESSIONAL USE

Read completely before use:

Clarident X is a complete in-office teeth whitening treatment kit for 1 patient. Hydrogen peroxide is the whitening agent and the strongest oxidant among peroxides, producing cleaning and whitening of vital and non-vital teeth.

The effectiveness of Clarident is widely proven in stains caused by fluorosis, tetracycline, nicotine, some foods, age. However, in case of severe or third degree fluorosis and tetracycline stains, it is recommended to use the micro-abrasion technique with Clarident TA.

Instructions for use:

Make an accurate diagnosis of the origin of the abnormal coloration.

It is not necessary to anesthetize.

1. Check that the gum is healthy in the area to be bleached.
2. Treat caries and possible leaks.
3. Treat sensitivity before applying the bleach.
4. Remove bacterial plaque. It is advisable to carry out a complete treatment prophylaxis. Do not use materials that contain glycerin or fluorine.
5. Isolate with rubber dam the teeth to be whiten or protect gums with gingival barrier.
6. Evaluate the shade of the tooth before treatment and, if possible, take a photo.
7. Prepare the whitening gel as follows: Uncap the Factor B syringe and fit the connector tightly. Then uncap the Factor A syringe and attach it to the connector. Pass Factor A into the syringe with Factor B and vice versa. The passage is repeated 6 to 8 times until homogenized. Finally remove the connector and the empty syringe and attach the applicator tip.
8. Apply the gel with the applicator tip on the teeth to be bleached with a thickness of 2 to 3 mm. Do not use metal spatulas.
9. Let it act for 40 minutes.
10. Suck the gel, rinse with plenty of water, remove the gingival isolation and finally the patient will have to rinse the mouth.
11. Evaluate the tonality and record with photography.
12. The number of sessions will depend on the severity of the case. At the end of the sessions make a topicalization with KDESIN, a desensitizer based on potassium nitrate and fluorine. If the patient has sensitivity, it will disappear in one or two days.

Recommendations:

- Wear gloves
- During the treatment, the enamel is especially permeable.

Warn the patient about the consumption of cigarettes or foods that can

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pigment it. If the treatment was carried out in one session, the white diet should be carried out for at least two days. If the product was stored in a refrigerator, use it after you have purchased room temperatura.

Caution:

Do not put in contact with the skin or mucosa. If this occurs, wash with plenty of water for a few minutes and cover the lesion with TEDEQUIM Gingival Protector.

Contraindications:

- Do not use in patients with dentin sensitivity.
- Do not apply to dentin or decayed teeth.
- Do not expose in areas with possible leaks.

Presentation:

- Factor A: hydrogen peroxide, 1.5 ml
- Factor B: gel 0.5 g
- KDESIN Desensitizing agent based on potassium nitrate and sodium fluoride, 2 ml
- Accessories.

Storage:

Keep between 4 and 15 C.

Technical Direction: Miryam Della Vedova - Farmacéutica - Mat. 2107

Authorized by ANMAT PM 1277-29

Exclusive sales professionals and healthcare institutions.

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Rev. A

KDESIN Dentin Desensitizer

KDESIN is a desensitizing agent containing potassium nitrate and sodium fluoride for treatment of dentinal hypersensitivity. The efficiency of potassium nitrate as a desensitizing agent is well known for the long time use of tooth's paste, but a better and more complete treatment can be achieved by combining both the potassium's and fluoride's effects.

KDESIN contains a high concentration of potassium nitrate, which offers a quick alternative for a significant reduction of sensitivity in a few days. It is used for root sensitivity cases caused by gingival retraction, in cases of sensitivity produced by temperature or chemical changes, dental bleaching, etc.

KDESIN is likely to be used by the dentist in office or the patient at home, applying once a day after a careful teeth brushing.

KDESIN should be indicated in the treatment of dentin hypersensitivity caused by acidogenic diets, bad habits, poor brushing techniques and the use of dental whiteners. Dentinal exposure can be due to processes, both physical and chemical, that lead to the loss of either enamel and cement or gingival tissue. Causal factors rarely act in isolation and include erosion, attrition, abrasion, bruxism, whitening, medication, aging, genetic conditions, gingival recession, and periodontal disease.

Instructions for use:

1. Remove dental plaque. If patient uses KDESIN at home, encourage thorough brushing.
2. In the case of use after Dental Bleaching, place gel in the same tray that was used with bleaching. In other cases place with the plastic needle or brush.
3. Apply KDESIN Gel surrounding the sensitive area. Let KDESIN remain in place for a few minutes or as much as possible. Do not wash. Avoid meals for at least one hour.
4. Remove by suction.

Cautions, Warnings and Recommendations:

- Due to the fact that the product acts when the ions K^+ and F^- penetrate, it is convenient to let it act as long as possible.
- In case the patient uses KDESIN at home, it is better to apply it after brushing at night and not to rinse.
- Avoid meals for at least one hour after application.
- For use in children between 6 and 12 years old, consult with the dentist.
- Keep away from children.
- Do not ingest.

- If you are taking any medication or are pregnant or nursing, consult your doctor before using this product.
- If you have any doubts, consult your dentist.
- If it is used in whitening cases and the patient has sensitivity, it is convenient to treat the sensitivity with KDESIN before whitening.

Restrictions:

Do not use in children under 6 years old.

Composition

15% potassium nitrate, 1% sodium fluoride, 84% inert ingredients.

Presentation:

Gel: 2 ml syringe, 1 or 3 units

Storage:

Keep between 4° and 30° C.

Technical Direction: Miryam Della Vedova – Farmacéutica - Mat. 2107

Authorized by ANMAT PM 1277-18 Exclusive sales professionals and healthcare institutions.

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