**Tedequim**

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**TEDEQUIM® CLARIDENT AUTO Y FOTOACTIVABLE**

**PÁGINA..... 2; 3; 4; 5; 6**

**INSTRUCTIONS FOR USE**

**TEDEQUIM® TEETH WHITENING KIT**

**PAGE ..... 7; 8; 9; 10; 11**

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**[www.tedequim.com](http://www.tedequim.com)**

**Rev. C**

**Clarident**  
**Teeth Whitening Kit**  
**Internal and External**  
**Hydrogen peroxide 35 or 38%.**  
**Self or photo activated**  
**Read completely before use.**  
**In-office use**

**Clarident** is a kit for complete in-office teeth bleaching treatment.

Clarident's bleaching agent is hydrogen peroxide, being the strongest oxidant among peroxides, producing cleaning and bleaching of vital and non-vital teeth.

The Oxigel powders that are included in the accessories contain sodium perborate and other components that make it possible to use the products with or without light. For photoactivated bleaching, any bleaching lamp emitting between 420 and 530 nm can be used.

The Kits are designed for external or internal use.

The effectiveness of Clarident is widely proven in stains produced by fluorosis, tetracycline, nicotine, some foods, age and endodontic trauma. However, in the case of severe or third degree fluorosis and tetracycline stains, it is advisable to use the microabrasion technique with **Clarident TA**.

**Instructions for Use:**

Make an accurate diagnosis of the origin of the abnormal coloration.

It is not necessary to anesthetize

**VITAL TEETH: EXTERNAL USE**

1. Control that gums are healthy in the area to be whitened.
2. Treat caries and possible filtrations.
3. Treat sensibility before the bleaching application.

**Bleaching lamp treatment:**

4. Remove bacterial plaque. It is convenient to perform a complete prophylaxis treatment. Do not use materials containing glycerin or fluoride.

5. Cover the soft tissue (gums and lips) with GINGIVAL PROTECTOR.

6. Isolate with rubber dam the teeth to whiten or protect gums with TEDEQUIMMR Gingival Barrier.

7. Cover again with GINGIVAL PROTECTOR the tissue that may be exposed after the previous procedure. It should be considered that hydrogen peroxide in these concentrations is very irritating for the mucosa.

8. Evaluate the shade of the tooth before treatment and if it is possible take a photograph.

9. Place a safety glasses on the patient.

10. Etch with acid (phosphoric up to 37%) the surface to be bleached for 10 seconds (only in cases of deep and specific stains). Then wash thoroughly and dry. This operation is not strictly necessary, it can be done only in the first session or it can be avoided. It is at the discretion of the professional.
11. Prepare the bleaching gel as indicated below.
12. Apply the gel with the help of a spatula or brush on the teeth to be whitened with a thickness of 2 to 3 mm. Do not use metallic spatulas.
13. Place the lamp close to the mouth and choose the intensity and time according to the professional's criteria. The total time of illumination is approximately 40 minutes distributed in two cycles of 20 minutes respectively or what is indicated by the manufacturer of the lamp.
14. When the illumination time is completed, remove the lamp and suction the gel and if it is necessary finish removing it with cotton.
15. Repeat items 11 to 14 until the wished result is obtained.
16. Suction the gel, rinse with plenty of water, remove the gingival isolation and ask the patient to rinse his mouth.
17. Evaluate the shade and record with a photograph.
18. The number of sessions will depend on the severity of the case. At the end of the sessions perform a topicalation with fluoride at neutral pH, Neutrodent or with KDESIN based on potassium nitrate and fluoride. If the patient has sensitivity this will disappear in one or two days.

#### **Non-bleaching lamp treatment:**

4. Remove bacterial plaque. It is convenient to perform a complete prophylaxis treatment. Do not use materials containing glycerin or fluoride.
5. Cover the soft tissue (gums and lips) with GINGIVAL PROTECTOR.
6. Isolate the teeth to be whiten with rubber dam or protect gums with Gingival Barrier TEDEQUIM-MR. 7- Cover again with GINGIVAL PROTECTOR the tissue that may be exposed after the previous procedure.
7. Cover again with GINGIVAL PROTECTOR the tissue that may be exposed after the previous procedure. It is important to consider that hydrogen peroxide in these concentrations is highly irritating to the mucosa.
8. Evaluate the shade of the tooth before the treatment and take a photograph if possible.
9. Etch with acid (phosphoric up to 37%) the surface to be bleached for 10 seconds (only in cases of deep and punctual stains). Then wash and dry thoroughly. This operation is not strictly necessary, it can be done only in the first session or can be avoided. It is at the discretion of the professional.
10. Prepare the bleaching gel as indicated below.
11. Apply the bleaching gel with a spatula or brush on the teeth to be whitened with a thickness of 2 to 3 mm. Do not use metallic spatulas.
12. Leave the bleaching gel working until it discolors or dries in about 20 minutes.
13. Suction the gel and if necessary finish removing it with cotton.
14. To finish the bleaching in the same session you can repeat items 10 to 13 or

place a few drops of peroxide on the same gel to moisten it again and continue.  
15. Suction the gel, rinse with plenty of water, remove the gingival isolation and then the patient should rinse his mouth.

16. Evaluate the shade and record with a photograph.

17. The number of sessions will depend on the severity of the case. At the end of the sessions perform a topicalization with fluoride at neutral pH, Neutrodent or with KDESIN based on potassium nitrate and fluoride. If the patient has sensitivity this will disappear in one or two days.

#### **NON-VITAL TEETH: INTERNAL USE**

Internal bleaching is an effective way to lighten teeth stained by endodontic trauma, however it is limited to non-vital teeth.

This product should not be used on teeth with silver cones. It is important to first remove the silver cones, treat the roots with gutta percha and wait at least one month before starting treatment with Clarident.

#### **Treatment:**

1. Check that the gingiva is healthy in the area to be bleached.
2. Evaluate the shade of the tooth before the treatment and if possible take a photograph.
3. Cover the soft tissue with GINGIVAL PROTECTOR (gingiva and lips).
4. Isolate the tooth to be treated with rubber dam.
5. Cover again with GINGIVAL PROTECTOR the tissue that may be exposed after the previous maneuver. It should be considered that hydrogen peroxide in this concentration is highly irritating to the mucosa.
6. Remove the filling material up to the base of the pulp chamber and seal.
7. Perform a complete prophylaxis treatment, without using materials containing glycerin or fluoride.
8. Etch the walls of the chamber with phosphoric acid and then rinse well.
9. Prepare the bleaching gel as indicated below.
10. Place the bleaching gel inside the pulp chamber and cover the external sides as well.
11. The action of Clarident is accelerated by heat or light. Light is more recommended. Light for about 15 minutes or until the product dries and discolors. In the absence of light, leave it act chemically until the product discolors or dries.
12. Rinse well and place the gel (or a cotton swab soaked in the peroxide) back inside the pulp chamber, provisionally obturate and leave it until the next visit. With time the concentration of peroxide will decrease and in addition to acting as a bleaching agent, it disinfects.
13. Items 9 and 10 can be repeated in the same session.
14. Evaluate the shade and record with a photograph.
15. The number of sessions will depend on the severity of the case. At the end of the sessions perform a topicalization with fluoride at neutral pH, Neutrodent or with KDESIN based on potassium nitrate and fluoride. If the patient has sensitivity this will disappear in one or two days.

16.To restore with composite or ceramic it is recommended to wait 1 to 6 weeks after bleaching to achieve a good bond strength and adhesion.

#### **Preparation of the bleaching gel**

Place the Oxigel powder in a Dappen cup (approx. one measure per element) then add drop by drop the hydrogen peroxide spatulating until a homogeneous gel is formed.

Use immediately as peroxide decomposes on exposure to air, light or heat.

\*NOTE: The dosage to be prepared will depend on the number of elements to be treated.

Use immediately.

It is recommended to close the container immediately after pouring the liquid.

#### **Recomendations:**

- Use gloves.
- During the treatment the tooth enamel is especially permeable. The patient should be warned about the smoking or eating foods that can pigment the enamel. If the treatment was performed in one session, the white diet should be done for at least 2 days.
- Store the bottles of hydrogen peroxide in the refrigerator and open one at a time. Close them immediately. Use the peroxide after it has reached room temperature.

#### **Caution:**

Do avoid contact with skin or mucous membranes, in case of contact, wash with abundant water for 15 minutes and cover the lesion with Gingival Protector.

#### **Contraindications:**

Do not use in patients with dentin sensitivity.

Do not apply on dentin or carious teeth.

Do not expose in areas with possible filtrations.

#### **Presentation:**

Clarident PH 35%

Clarident PH 38 % photoactivated

#### **The kit contains:**

- Gingival and lip protectant: 9 g
- Hydrogen peroxide 35% or 38%: 10 ml
- Oxigel powder: 2.5 g
- Accessories.

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