Tedequim

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Root Canal Chemical Stretcher

EDTA disodium 17%

EDTA stands for Ethylenediaminetetraacetic acid. It is a hexadentate chelating agent. The main chemical property of EDTA and the one that justifies its use in dentistry is to trap cations, such as Ca ++, Mg ++, Fe ++, Cu ++, Hg ++, Zn ++, etc. forming a soluble complex. Actually, the dissociated acid is the one that intervenes in the chelation reaction. Generally the disodium salt of EDTA is used at physiological pH ie at pH = 7.4.

Among the chelating agents EDTA is the most important to trap calcium ions, therefore EDTA has been used for a long time and is currently applied in dental therapeutics in those situations in which it is necessary to remove dentine or smear layer, either in conservative dentistry or in endodontics. So in the root canal biomechanical preparation, when it is desired to achieve chemical widening in an innocuous and simple way and to facilitate the expansion of narrow canal it is advisable to use EDTA solutions.

In order to gather the best properties of this solution, it has been proposed to associate disodium EDTA with sodium hydroxide to achieve the optimum pH and a cationic quaternary ammonium surfactant agent. This is achieved by increasing the bactericidal power of the product, improving the EDTA diffusion, lowering the surface tension and increasing the cleaning power on the wall by reducing the smear layer or viscous layer adhered and accumulated during the instrumentation.

The advantages offered by this product in the preparation of the root canal are the following:

- Chemical stretching by superficial softening of the dentin wall.
- Elimination of the residual layer.
- Better mechanical cleaning of the dentin wall.
- Disinfection or antibacterial action of the dentin wall.
- Increased dentinal permeability to medications.
- Greater adhesion of the cement to the dentin wall.

Instructions for use:

The application must be done by placing a few ml. in a Dappen glass, wetting the instruments and inserting it by pumping it into the root canal. It can be repeated as many times as necessary. It can also be placed with a syringe irrigating the canal. A constant renewal of the product is more efficient than a continuous application.

Form of presentation:

- Dropper bottle x 25 ml.
- Glass bottle x 20 ml.

Composition:

EDTA Na ₂ . 2 H ₂ O	18.8 g.
Quaternary ammonium cationic detergent	1.0 g.
Sodium hydroxide: csp reach	pH 7.4
Bidistilled water	csp 100 ml.

Contraindications:

They have not been found so far.

Precautions:

Use with caution in case of stepladders or apical blocks where instruments are forced since by softening the dentin wall a perforation could occur.

If leakage through the apical foramen towards the tissues of the peri-apex occurs during the preparation of the root canal, a decalcifying and irritative action of the periapical bone can be exerted that will heal in a few days without affecting non-calcified tissues. It is convenient to wash after application.

Recommendations:

When used as irrigating agent it is recommended to alternate with sodium hypochlorite solution.

Storage:

Keep the container tightly closed. Store at temperatures between 4° and 40° C.

Approved by ANMAT PM 1277-12

Technical direction: Miryam Della Vedova - Pharmacist - Licence 2107

Exclusive sale to professionals and sanitary facilities

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