Tedequim

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Rev. A

TEDEQUIM® GINGIVAL BARRIER

The gingival barrier is a photocurable methacrylate-based resin for gum protection in Professional Teeth Whitening or to isolate the tissue adjacent to the tooth to be treated in endodontics.

Provides good adhesion and easy removal of it at the end of the treatment. It hardens with either halogen light curing lamps, or with LEDs between 410 and 500 nm. It stands out for having a low release of heat when polymerizing and for having little bitter taste that is well tolerated by the patient.

Instructions for use:

For vital whitening

- 1) Remove bacterial plaque.
- 2) Thoroughly dry teeth and gums.
- 3) Uncap the syringe and place the applicator needle.
- 4) Apply the barrier about 4 to 6 mm wide on the gum. Seal the interproximal spaces. Place the resin approximately 0.5 mm on the dry enamel to seal.
- 5) Apply the resin on a tooth beyond the last one to be bleached.
- 6) Make sure that all the fabric has been covered.
- 7) Light cure the resin for 20 seconds for each guide width.

Non-vital whitening

- 1) Remove bacterial plaque.
- 2) Thoroughly dry teeth and gums.
- 3) Uncap the syringe and place the applicator needle.
- 4) Apply the barrier about 4 to 6 mm wide on the gingiva, lingual and labial, and between the interdental areas.
- 5) To avoid the heat of the curing lamp, cover approx. the adjacent tooth halves with a 3-4mm layer of barrier.
- 6) Light cure the resin for 20 seconds for each width of guide to polymerize.

For endodontic treatments

- 1) Remove bacterial plaque.
- 2) Thoroughly dry teeth and gums.
- 3) Place Rubber Dam.
- 4) Uncap the gingival barrier syringe and attach the applicator needle.
- 5) Fill tissues that may be exposed with gingival barrier.
- 6) Light cure the resin for 20 seconds

Gingival Barrier Removal

The Light Cure Gingival Barrier is quickly and easily removed in one or more pieces. Check the interproximal areas so that no residue remains. If whitish

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areas appear (damaged by peroxide), apply Clarident Gingival and lip protector to heal the lesion.

Collateral or secondary actions:

They have not been found so far.

Precautions:

- Use a mouth opener
- Follow the manufacturer's instructions for the curing light you use.

Warnings and Restrictions:

- Do not apply to patients with ulcerative gingivitis
- Do not apply to patients allergic to resins.
- Do not eat.
- Keep away from children.

Cautions in handling:

Cap immediately after using the syringe. Prevent the curing light from reaching the syringe. Keep in a cool place.

The product is sensitive to light therefore it must be kept well covered and inside the box to protect from light.

Presentation: 1.2 ml syringe

Storage:

Store well covered, protected from light and at temperatures between 15 and 30°C preferably.

Autorized by ANMAT PM 1277-25

Technical director: Miryam della Vedova - Pharmacist - Licence 2107

Exclusive sale to professionals and sanitary facilities

**Manufactured by TEDEQUIM SRL
Bv. De los Polacos 6136- Los Boulevares
X5147GGP Córdoba - Argentina
Tel/FAX: 54-3543-448260
www.tedequim.com.ar**