Tedequim

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www.tedequim.com.ar

Rev. C


Clarident PC
Dental bleaching. Home use
Urea peroxide 16% - 22%

Dental whitening is a simple and conservative treatment within aesthetic dentistry, which can return the tooth its color and translucency, when it presents pigmentations produced by extrinsic or intrinsic factors.

It is a non-invasive treatment that must be used by the professional without risk because it:

- Preserves the tooth structure
- Improves the aesthetics of the patient
- Simplifies later restorative procedures.
- Improves oral health because patients who receive cosmetic procedures such as teeth whitening have an improvement in oral hygiene habits and general oral health.

There are different whitening products in different presentations and concentrations. However, the active agent is the same: hydrogen peroxide. The carbamide peroxide or urea peroxide is a source of hydrogen peroxide since it is a 1:1 complex between hydrogen peroxide and urea:

Urea peroxide

Hydrogen peroxide + urea

This product is a gel that must be applied on the enamel by means of a tray, made by the professional, which is adapted to the mouth of the patient.

The result of the treatment and the number of times it should be applied will depend on the good preparation of the tray, since the better the contact of the gel with the enamel, the better the performance.

This ambulatory technique is performed in patients with constant behavior who have not very pronounced teeth color changes or as a complement to the procedure in the office. It can also be used in maintenance cases.

Instructions for use:

1. Make an accurate diagnosis of the origin of the abnormal coloration. It is not necessary to anesthetize
2. Check that the gum is healthy in the area to be bleached.
3. Treat caries and possible leaks.
4. Treat sensitivity problems before treatment. (KDESIN-Neutrodent)
5. Take a picture before treatment and use a color chart.
6. Remove bacterial plaque. It is convenient to carry out a complete prophylaxis treatment. Do not use materials that contain glycerin or fluoride.
7. Making the tray: Make a model of the hemiarcade in the shape of a horseshoe. A photocurable spacer resin is placed on the obtained model. The

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cuvette is made in the clinic or in the laboratory with a vacuum and heat machine (vacupress or similar). This will be done in soft plastic from 0.6 to 0.9 mm thick. It should cover all the teeth of the arch and extend 2 to 3 mm towards the apical gingival margin. In the upper arch, the palate should be left free and spread over the palatal wrinkles.

8. Remove the syringe cap and place the plastic applicator. There is no need to remove the applicator after use, the cap can also be placed in its tip.

9. Introduce the gel in the spaces of the tray corresponding to the vestibular surface of each tooth to be bleached. The patient should be prevented from placing excess gel because it can irritate the gum.

10. Place the tray or cuvette with gel in the corresponding arch and leave it for at least two hours.

11. Remove the cuvette, sanitize it with a toothbrush and make a careful oral cleaning.

12. Repeat steps 7 to 10 as many times as necessary until the color remains stable.

13. 48 hours after the treatment is finished, the pieces adopt the final color.

14. If it is desired to restore with composite or ceramic, it is advisable to wait 1 to 6 weeks after whitening, to achieve a good bonding strength and adhesion.

Recommendations:

- Use gloves.
- During the treatment, the enamel is especially permeable. The patient should be warned about the consumption of cigarettes or foods that can pigment it. To prevent this, perform a '*white diet*' during the entire treatment and 4 days after the it is concluded.
- Store the syringes in the refrigerator and open one at a time. Close it immediately after use.
- Keep out of the reach of children.
- In the case that the patient has sensitivity, use KDESIN or Neutrodent.

Cautions:

- Prevent contact with skin or mucosa. In case this happens, wash with plenty of water.
- If there are restorations in anterior teeth, they should be replaced once the treatment is finished to match the color.
- If there are fixed porcelain prostheses in anterior teeth, they must be replaced once the treatment is finished to match the color.

Contraindications:

You should not bleach if there are:

- Gingival pathologies, periodontal diseases, such as gingival retraction. In this last case it is highly advisable to do a good cleaning and reduce inflammation in gums first.

- Dentin hypersensitivity. It is convenient to treat the sensitivity before and then do the whitening.
- Poor restorations.
- Caries and possible leaks.
- Young teeth.
- During pregnancy.
- Teeth with deep fissures.
- Teeth with fine enamel thicknesses - Significant loss of enamel.
- Dental trauma.
- Radicular reabsorption.
- Enamel development defects.
- Pigmentation caused by corrosion of amalgams.
- Composites poorly adjusted.
- Teeth with large restorations or repeated seals in the same piece.
- Very dark teeth.
- Abnormal dental morphology (its internal structure may be rare).
- Tetracycline grade III and moderate and severe fluorosis.

Restorations and porcelain pieces will not whiten.

Its white fillings and crowns will not bleach but will get clean.

This can cause the teeth to be noticeably whiter than the reconstructions.

Most patients prefer to change the dark reconstructions to make the rest of the teeth whiter and brighter. This will ensure that the best tonality is achieved. It is recommended to perform the whitening before treating the teeth, waiting 20 days.

Presentation:

Single syringe x 3 g

Storage:

Store preferably in a refrigerator. Use at room temperature.

This product must be used exclusively under the supervision of the professional

Authorized by ANMAT PM 1277-29

Technical director: Miryam Della Vedova - Pharmacist - License 2107

Manufactured by TEDEQUIM SRL

Bv. De los Polacos 6136, Los Boulevares.

X5147GGP Córdoba, Argentina.

Tel/FAX: 54-3543-448260

www.tedequim.com

Rev. C