Tedequim

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TEDEQUIM®
EDTA EN GEL 24%

EDTA is the abbreviation for ethylenediaminetetraacetic acid. It is a hexadentate chelating agent. The main chemical property of EDTA and the one that justifies its use in odontology is to trap cations, such as Ca⁺⁺, Mg⁺⁺, Fe⁺⁺, Cu⁺⁺, Hg⁺⁺, Zn⁺⁺, etc. forming a soluble complex. Among the chelating agents EDTA is the most important for trapping calcium ions, therefore EDTA has been applied for a long time and is currently applied in dental therapeutics in those situations in which it is necessary to remove dentin or dentinal mud, either in periodontics or endodontics. The gel is water soluble allowing its complete elimination with water or irrigation solutions in aqueous medium. The concentration of 24% does not promote inconvenience for dentin or periapice, but confers to the product a faster and more efficient action. It takes 2 to 3 minutes until the product take effect. At physiological pH the predominant species is trisodium EDTA.

Instructions for use:

1. Allow the product to come to room temperature before use.
2. Remove the cap of the syringe and insert the needle applicator.

a. Endodontics:

- 3a. Place a small amount of EDTA in the entrance of the root canal and begin instrumentation.
- 4a. Irrigate and if necessary repeat the operation by placing new amounts of the product.
- 5a. After finishing the instrumentation irrigate the canal with Sodium Hypochlorite or the irrigating solution considered convenient before obturation.

b. Periodontics:

- 3b. Remove granulation tissue.
- 4b. Apply EDTA gel on the exposed root surface. Leave for 2 to 3 minutes and perform root curettage.
- 5b. Repeat the operation if necessary.
- 6b. Rinse with physiological solution or demineralized water.
- 7b. Suture if necessary.

Presentation:

- 1 syringe x 3 g and 3 applicator tips.
- Kit containing 3 syringes x 3 g and 9 applicator tips.

Composition

Ethylenediaminetetraacetic acid disodium salt dihydrate, sodium hydroxide, distilled water and thickener.

Contraindications:

No contraindications have been found so far, however it should not be used in people who report sensitivity to any component of the product.

Side effects:

The product does not cause necrosis in surrounding periodontal tissues.

Precautions, warnings and special cares:

- Use with caution in case of apical steps or blockages, where the instruments are forced, since by softening the dentin wall, perforation could occur.
- If during the preparation of the canals there is leakage through the apical foramen into the periapical tissues, a decalcifying and irritating action of the periapical bone can be exerted, which will heal in a few days without affecting non-calcified tissues. It is advisable to wash after application.
- Do not use the product after the shelf life period.
- Do not use if the container is damaged.
- Dispose of the product according to the local legislation in force.
- Discard the applicator needle after use.
- Keep the product out of the reach of children.

Storage:

Keep the product between 4 and 30 °C. Use at room temperature.

Technical Direction: Miryam Della Vedova - Pharmacist - Mat. 2107

Authorized by ANMAT PM 1277-33

Exclusive sale to professionals and health institutions.

Argentine Industry

Elaborated by:

TEDEQUIM S.R.L.

Bv. De los Polacos 6136, Los Boulevares.

X5147GGP, Córdoba, Argentina

TEL: 0054-3543-48260

www.tedequim.com

Rev. A