Tedequim

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FAgamin®

PÁGINA..... 2; 3

INSTRUCTIONS FOR USE

FAgamin®

PAGE..... 4; 5

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Bv. de los Polacos 6136, X5147GGP


Córdoba, Argentina

FAgamin®
Antibacterial, remineralization, desensitizing
agent to inhibit and arrest dental caries.

Professional use only.

Description:

FAgamin® contains Silver diamine fluoride (Ag (NH₃)₂F). Other similar products used in caries prevention and dental hypersensitivity have several disadvantages in opposition to **FAgamin®**. When applying silver nitrate (AgNO₃) or sodium fluoride (NaF), the dental enamel hydroxiapatite dissociates resulting in soluble calcium and phosphate ions. The action on teeth of silver diamine fluoride is based on the chemical reaction with hydroxiapatite.

Calcium fluoride (CaF₂) and silver phosphate (Ag₃ PO₄) form deposit over the tooth surface, reinforcing and protecting it. This explains the remineralizable effect of **FAgamin®**.

The silver ion (Ag⁺) inhibits the development of bacteria and has an anti-enzymatic action. The ammonia stabilizes the silver ions by forming a complex and keeping them in aqueous solution, giving a basic pH. On the other hand, the fluoride ion decreases the solubility of hydroxyapatite and inhibits the development of already formed caries. In case of applying **FAgamin®** on a cavity, a black deposit is formed due to the reaction of silver ions with organic matter. That's why it also acts as a caries detector.

FAgamin® has a wide range of clinical applications due to its chemical composition. Some of these are: old teeth enamel and caries nearby dentin remineralization, caries detection and arrest, incipient and old caries inhibition, dental hypersensitivity desensitization. It also has bactericidal and antiseptic action.

Furthermore, **FAgamin®** offers a practical, easy and quick topical method for dental caries treatment on temporary teeth.

Applications:

It is applied in temporary teeth, in incipient or rampant caries, in pits and fissures (after widening) if the dentist deems it convenient, in cavity preparations for prosthesis, crowns, dies, fixed bridges. It can be applied in root canals diluting it (1:10 with distilled water).

In patients with severe cognitive and physical disabilities and patients with dental phobias.

Instruction for use:

Wash and dry the tooth before starting treatment.

FAgamin® is applied using a swab or brush impregnated with 2 or 3 drops of the product. It is placed on the area of the tooth to be treated, left to act for a few minutes by rubbing the applied area. Finally, rinse well to remove traces of the solution.

To avoid the contact of **FAgamin®** with the mucous membrane, isolate it with a cotton pellet or rubber dump.

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In cases where the carious tissue is near gingival tissue, cover it with vaseline or cacao butter.

In incipient caries, **FAgamin®** is applied only once, but in other cases more than one application may be necessary.

Handling care:

- Attach the bottle cap immediately after using.
- Wear gloves.
- Store in a place protected from light.
- In case of staining the skin, wash it with household bleach rubbing with cotton, or hydrogen peroxide.

Warnings:

- Due to the fact that the reaction of the silver on the dentin can produce darkening, it is not advisable to use **FAgamin®** in permanent front teeth where the aesthetic factor is important.
- To avoid the aesthetic factor, the treated area can be covered with glass ionomer or potassium iodide. The immediate application of a potassium iodide solution reduces the darkening produced by the reaction of silver ions.
- In case of deep cavities, since filtrations to the pulp may occur, it is advisable to avoid its use or dilute the dose with distilled water, in a 1:10 ratio.
- Leakage to the pulp may also occur depending on the quality of the tooth.
- When contact occurs accidentally or when the patient manifests pain, rinse immediately with mouthfuls of salt water or water with 10 V hydrogen peroxide. If pain persists, apply camphorated phenol for dental use.

Contraindications:

Generally does not possess in application doses. Do not apply too close to pulp or when there are known deep seepages.

Content:

FAgamin® contains 5 ml of 38% silver diamine fluoride solution at controlled pH

Storage:

Keep at temperatures between 4 and 30 °C, well closed and protected from light.

Exclusive Professional use.

Technical director: Miryam della Vedova – Pharmacist – Licence 2107

Autorized by ANMAT PM 1277-14

Made in Argentina

Manufactured by TEDEQUIM SRL

Bv. De los Polacos 6136- Los Boulevares

X5147GGP Córdoba – Argentina

Tel/FAX: 54-3543-448260

www.tedequim.com.ar

Rev. D