**Tedequim**

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## **KDESIN**

### **Dentin desensitizer**

KDESIN is a desensitizing agent containing potassium nitrate and sodium fluoride for treatment of dentinal hypersensitivity. The efficiency of potassium nitrate as a desensitizing agent is well known for the long time use of tooth's paste, but a better and more complete treatment can be achieved by combining both the potassium's and fluoride's effects.

KDESIN contains a high concentration of potassium nitrate, which offers a quick alternative for a significant reduction of sensitivity in a few days.

It is used for root sensitivity cases caused by gingival retraction, in cases of sensitivity produced by temperature or chemical changes, dental bleaching, etc.

KDESIN is likely to be used by the dentist in office or the patient at home, applying once a day after a careful teeth brushing.

KDESIN should be indicated in the treatment of dentin hypersensitivity caused by acidogenic diets, bad habits, poor brushing techniques and the use of dental whiteners. Dentinal exposure can be due to processes, both physical and chemical, that lead to the loss of either enamel and cement or gingival tissue. Causal factors rarely act in isolation and include erosion, attrition, abrasion, bruxism, whitening, medication, aging, genetic conditions, gingival recession, and periodontal disease.

#### **Instructions for use:**

##### **KDESIN Liquid:**

- 1) Remove bacterial plaque. If the patient uses it at home, make a good brushing.
- 2) Soak a cotton swab, swab or brush with 2 or 3 drops of the liquid with KDESIN
- 3) Apply on the sensitive parts of enamel, dentin or cement.
- 4) Let it act as long as possible. It is advisable not to rinse and not eat for at least an hour.

##### **KDESIN Gel:**

- 1) Remove bacterial plaque. If the patient uses it at home, make a good brushing.
- 2) In the case of teeth whitening, place with the help of the applicator a little gel in a tray that is used for whitening. In other cases it can be applied with a brush or brush.
- 3) Apply on the sensitive parts of enamel, dentin or cement.
- 4) Leave to act for a few minutes. Do not eat for at least an hour.
- 5) Remove by suction or rinse.

### Cautions, Warnings and Recommendations:

- Due to the fact that the product acts when the ions K<sup>+</sup> and F<sup>-</sup> penetrate, it is convenient to let it act as long as possible.
- In case the patient uses KDESIN at home, it is better to apply it after brushing at night and not to rinse.
- Avoid meals for at least one hour after application.
- For use in children between 6 and 12 years old, consult with the dentist.
- Keep away from children.
- Do not ingest.
- If you are taking any medication or are pregnant or nursing, consult your doctor before using this product.
- If in doubt, consult your dentist and/or pharmacist. If it is used in whitening cases and the patient has sensitivity, it is convenient to treat the sensitivity with KDESIN before whitening.

### Restrictions:

Do not use in children under 6 years old.

### Handling care:

The product is an aqueous solution of stable inorganic salts that does not require special handling care.

It is only advisable to cover the container after handling.

### Composition:

15% potassium nitrate, 1% sodium fluoride, 84% inert ingredients.

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