**Tedequim**

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**Rev. B**

**TEDEQUIM®**  
**ASTRINGENT SOLUTION**  
**Without epinephrine**

It is possible to retract the gingival tissue without risk using cotton cords soaked in a hemostatic solution of ferric sulfate free of epinephrine, preventing the possible side effects that this constrict vessel could probably produce. The result lesion of the tissue is insignificant and it is likely to immediately repeat the cord packing if necessary. This packet pushes laterally the gingival tissue and exposes the dentinal surface between 0.35 and 0.5 mm. beyond the gingival line of the sulcus to assure an accurate impression or crown restorations.

This solution can be used with the cotton cords of the TEDEQUIM kit, or with any other cord.

**Instructions of use:**

The cotton cord size is chosen depending on the sulcus depth.

The cord is cut and soaked in the ferric sulfate hemostatic solution using the dropper. DO NOT submerge the cord directly into the bottle. Then, the tooth's neck is rounded with the cord and packed, with an appropriate instrument, where the sulcus is under the level of the gingival crest. When the cord is properly set, it is dampened again using a swab.

After 3 or 5 minutes, the sulcus is observed in order to check the retraction. If the result was not the one expected, another cord packing is likely to be done with new material or added astringent solution.

When removing the packet, astringent solution is needed to wet the cord in order to prevent bleeding.

The Astringent solution can be used any time when hemorrhage occurs. A soaked swab is pressed over the bleeding area until flow ceases.

This solution can be use with any cotton cord.

**Restrictions:**

Do not apply in patients with:

- Allergy.
- Ulcerative gingivitis.
- Stomatitis
- Acute hemorrhages

**Adverse reaction:**

No adverse reactions detected.

**Composition:**

15.5% ferric sulfate ( $\text{Fe}_2(\text{SO}_4)_3$ ) solution.

**Presentation:**

20 ml astringent solution

**Storage:**

Keep between 4°C and 40 °C protected from light.

**Autorized by ANMAT PM 1277-16**

**Technical director: Miryam della Vedova - Pharmacist - Licence 2107**

**Exclusive sale to professionals and sanitary facilities**

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