Tedequim

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www.tedequim.com.ar

Rev. B

TEDEQUIM^{MR}
ETCHING ACID FOR PORCELAIN
9, 5 - 10 % hydrofluoric acid buffered gel

Etching porcelain enhances the bond with composite, resins, cements and other organic resins. The gel etches the porcelain in 2 minutes. The product is indicated for etching porcelain veneers and inlays, fractured porcelain surfaces prior to bonding in bridges, crowns, etc.

Instructions of use:

1. Use protective eyewear and gloves.
2. Place rubber dam to isolate completely the operative site to protect soft tissue.
3. Remove syringe cap and place the plastic needle. Carefully apply the gel only in the area to be etched. An assistant must hold the vacuum system nearby the porcelain surface to evacuate the acid fumes.
4. Leave the gel for 2 minutes.
5. Remove the gel and rinse with copious amount of water.
6. Dry with compressed air.
7. Apply TEDEQUIM[®] Silane for better adhesion.

Cautions:

Hydrofluoric acid is a powerful and aggressive chemical and should be handled with care. In case of contact with soft tissue, skin or eyes, rinse the exposed area with large amount of water and neutralize with sodium carbonate acid or calcium gluconate immediately.

Use plastic or Teflon material for manipulation of this product. Do not use metal as it will corrode.

In case of direct contact with the skin or soft tissue, burns may be difficult to notice at the beginning. Fluoride ions penetrate rapidly causing soft tissue necrosis and decalcification of bones. In opposition to other acids, which are quickly neutralized, this process can continue for days.

Composition:

9.5- 10 % hydrofluoric acid in gel.

Presentation:

2 ml syringe and plastic needles

Storage

Keep between 4°C and 40 °C protected from light.

Autorized by ANMAT PM 1277-25
Technical director: Miryam della Vedova - Pharmacist - Licence 2107
Exclusive sale to professionals and sanitary facilities

Manufactured by TEDEQUIM SRL
Bv. De los Polacos 6136- Los Boulevares
X5147GGP Córdoba - Argentina
Tel/FAX: 54-3543-448260
www.tedequim.com.ar
Reb.B